

CERTIFICATE OF ANALYSIS

CODE: 38417

PRODUCT: VALQUER OXIDIZING CREAM 10 VOL (3%) 500 ML

AMOUNT: 2.995

BATCH: 240834

MANUFACTURING DATE: 07/06/2024

PHISICO-CHEMICAL ANALYSIS

PARAMETER	METHOD	SPECIFICATION	RESULT
TYPE OF EMULSION	IT-003-03	o/w	Pass
COLOR	IT-003-01	STANDARD	Pass
ODOR	IT-003-40	STANDARD	Pass
pH	IT-003-09	2 - 3	3

THE PRODUCT IS CONFORM TO ITS SPECIFICATIONS

DATE:

11/06/2024

SIGNATURE:



R.M. de Toledo nº 1, Tomo 1.513, Folio 180, Hoja T-32407

TECHNICAL DATA SHEET

Rev:

Revision Date:

Code:

38417

Description:

VALQUER OXIDIZING CREAM 10 VOL (3%) 500 ML

Development Number:

D01

Brand:

HAIR CARE

Ean 13:

8420212384170

Responsible person name:

VÁLQUER LABORATORIOS, S.L.

CPNP product reference:

3172855



LABELLING

TEXT

This oxidizing emulsion with hydrogen peroxide with a concentration up to 10 volume (3%) develops colour when mixed with colouring cream for permanent colourations. For hair bleaching, it is mixed with bleaching powder. Instructions for use: Mix the cream in a bowl with the colouring cream or the bleaching powder to form a homogenous cream. Apply on the hair.

VALQUER
LABORATORIOS

OXIDIZING CREAM

Oxidante en crema

10 VOL
(3%)

Professional

500 ml - 16.90 Fl. Oz.



REUSABLE
CONTAINER



TREATED
WATER



100%
VEGAN



DERMO
TESTED

VALQUER

LABORATORIOS

EN-OXIDISING CREAM 10 VOL (3%). Oxidising emulsion with hydrogen peroxide at a concentration of 10 volumes (3%) which develops the colour when mixed with the colouring cream in permanent colourings. For hair bleaching, it is mixed with bleaching powder. How to use: Mix in a plastic bowl with the colouring cream or bleaching powder until a homogeneous cream is formed and apply. Warning: Contains hydrogen peroxide. Use appropriate gloves. Avoid contact with eyes. Rinse eyes immediately if product comes into contact with them.

ES-OXIDANTE EN CREMA 10 VOL (3%). Emulsión oxidante con agua oxigenada a una concentración de 10 volúmenes (3%) que desarrolla el color al mezclarse con la crema colorante en las coloraciones permanentes. Para decoloraciones capilares, se mezcla con polvo decolorante. Modo de empleo: Mezclar en un bol de plástico con la crema colorante o polvo decolorante hasta formar una crema homogénea y aplicar. Advertencias: Contiene peróxido de hidrógeno. Utilizar los guantes apropiados. Evitar todo contacto con los ojos. Enjuagar inmediatamente los ojos si el producto entra en contacto con éstos.

FR- CRÈME OXYDANTE 10 VOL (3%). Émulsion oxydante avec de l'eau oxygénée à une concentration de 10 volumes (3 %) qui fait ressortir la couleur en la mélangeant à la crème colorante dans les colorations permanentes. Pour les décolorations capillaires, mélanger avec de la poudre décolorante. Mode d'emploi : Mélanger dans un bol en plastique avec la crème colorante ou la poudre décolorante jusqu'à l'obtention d'une crème homogène et appliquer. Mises en garde : Contient du peroxyde d'hydrogène. Porter les gants appropriés. Éviter tout contact avec les yeux. Rincer immédiatement les yeux si le produit entre en contact avec ceux-ci.

VPQ38417/05

PT-OXIDANTE EM CREME 10 VOL (3%)

DE-OXIDATIONS CREME 10 VOL (3%)

IT-OSSIDANTE IN CREMA 10 VOL (3%)

JP-パルカー酸化クリーム 10vol(3%)

CH- 氧化霜 10 体积 (3%)

AR- (3%) كريم فالكير المؤكسد مل

RU-КРЕМ-ОКСИД 10 VOL (3%).

FA- (3%) کرم اکسید کننده والکوتر 10

RO-CREMĂ OXIDANTĂ 10 VOL (3%)

PL-KREM UTLENIAJĄCY 10 VOL (3%)

BG-ОКСИДИРАЩ КРЕМ 10 VOL (3%)

INGREDIENTS: AQUA, CETEARYL ALCOHOL, HYDROGEN PEROXIDE, CETETH-24, SODIUM LAURYL SULFATE, PHOSPHORIC ACID, SODIUM STANNATE, TETRASODIUM PYROPHOSPHATE

FLACONE
PLASTICA
HDPE 2

TAPPO
PLASTICA
PP 5



12M



raccolla DIFFERENZIATA



VALQUER LABORATORIOS S.L.U.

P. I. La Dehesilla, 33. 45440 Villaminaya, Toledo, España.

VALQUER UK LTD

118 Greenwood Hill, Belfast, BT8 7WF, Northern Ireland, United Kingdom.

VALQUER LABORATORIOS MEXICO SA de CV

Calle Antiguo Camino a Xochimilco 35, Edificio D, Departamento 302

Colonia La Noria Xochimilco CP 16030. Ciudad de México, México.